

# PURELAB ACCESSORIES

## REMOTE DISPENSING AND MONITORING



# Bench-top Dispensers and Control Stations

Remote dispensing and monitoring enables you to position your PURELAB water purification system under the bench or on a wall to free up valuable space on the bench-top.

### Remote Dispense Gun

ELGA Dispense Guns offer fingertip control over the flow rate, of up to 2 liters per minute. Water purity is maintained at all times by recirculation through the Dispense Gun.

- Dispense Gun clips securely into a wall-mounting bracket for fixed dispensing or can be unclipped for extended reach (1 meter range from the bracket)
- The wall-mounting bracket can be positioned up to 2 meters from the water purification system

### Integral Dispense Gun

- Increases the dispensing range by 1 meter (3ft)

### Remote Control Station

Enables the user to monitor water quality and check system status from a convenient location without needing to make a visual check on the main system. The Remote Control replicates the display on the water purification system and can be positioned up to 3 meters away from the unit.

- Provides a choice of display modes for operational status via push button controls
- Wall mountable for flexibility and convenience
- RS232 port for connection to a printer for easy data collection for compliance with GLP guidelines



# PURELAB ACCESSORIES

## REMOTE DISPENSING AND MONITORING

### Specifications Remote Dispensers

Tube length	Approx 2 meters
Flow Rate(Depending on water system used)	0.8 – 2 l/min

### Ordering Information Remote Dispensers

Description	Part Number
Remote Dispense Gun – PURELAB Option-R 7/15 and Option-E 5/10	LA629
Remote Dispense Station (Including Dispense Gun) – PURELAB Ultra, Classic, Option-Q 7/15, Option-R 7/15 and Option-E 5/10	LA643
Integral Dispense Gun – PURELAB Ultra and Classic	LA644
Open Circuit Ring Main Kit – PURELAB Ultra and Classic	LA646
Closed Circuit Ring Main Kit – PURELAB Ultra and Classic	LA647

### Dimensions and Weight Remote Controls

Height	120mm
Width	200mm
Depth	70mm
Weight	1kg

### Ordering Information Remote Controls

Model	Part number
PURELAB Ultra and Classic	LA645
PURELAB Option-Q 7/15	LA675
PURELAB Option-S 7/15	LA624
PURELAB Option-R 7/15	LA625
PURELAB Option-E 5/10	LA626
PURELAB Prima 7/15/30	LA628



ELGA LabWater

Tel: +44 (0) 1494 887500 Fax: +44 (0) 1494 887505

Email: [info@elgalabwater.com](mailto:info@elgalabwater.com) Website: [www.elgalabwater.com](http://www.elgalabwater.com)

ELGA is the global laboratory water brand name of Veolia Water Solutions and Technologies. The information contained in this document is the property of VWS (UK) Ltd, trading as ELGA LabWater, and is supplied without liability for errors or omissions.

© VWS (UK) Ltd. 2010 – All rights reserved. ELGA® and PURELAB® are registered trademarks of VWS (UK) Ltd.

LITR38918-02